



# 22NRM01 TraMeXI

Traceability in Medical X-ray Imaging Dosimetry

*Summary of discussions on practical sessions*


*Workshop on X-ray imaging dosimetry  
20-22 November 2024 Helsinki University Hospital*

*Luigi Rinaldi  
INMRI-ENEA and Bambino Gesù Children's Hospital*



# Mammography

- ▶ **Measurement geometry:** Sechopoulos I, Dance DR, Boone JM, et al. Joint AAPM Task Group 282/EFOMP Working Group Report: Breast dosimetry for standard and contrast-enhanced mammography and breast tomosynthesis. Med. Phys. 2024;51:712–739. <https://doi.org/10.1002/mp.16842>



Distance correction: square law

$$K_{ref} = K_m \left( \frac{l_m}{500 \text{ mm}} \right)^2 C_k$$

Forward scatter correction factor

# Mammography

## ► Measurement:

- Ionization chamber:

Raw data (dosimeter readings) → calibration coefficient


→ temperature and pressure correction coefficient

- XMM:

Raw data (dosimeter readings) → calibration coefficient

- Can we compare them? **NO!**


IC reference point  $P_1$  and XMM reference point  $P_2$  are not at the same distance from the X-ray tube focus, use Inverse Square Law.



# Mammography

## ► Impact of compression paddle:

| Radiation Quality | In contact (0 cm) | Far (15 cm)     | Forward scatter factor |
|-------------------|-------------------|-----------------|------------------------|
|                   | Air kerma (mGy)   | Air kerma (mGy) |                        |
| W/Ag 28 kV        | 1.468             | 1.415           | 1.037                  |




When the paddle was adjusted to its maximum clinical distance, the proportion of the detected forward scattering was only 1 %.

The most consistent way of performing air kerma measurements is to position the compression paddle at the maximum distance from the dosimeter [1].



[1] Toroi P, Könönen N, Timonen M, Kortensniemi M. Aspects of forward scattering from compression paddle in the dosimetry of mammography. Radiat Prot Dosimetry 2013;154(4):439-45. <https://doi.org/10.1093/rpd/ncs257>.

# Mammography

## ► Positioning:




top view



## Interventional X-ray system

- ▶ Inherent filtration: 2.5 mm Al
- ▶ Field size: 10 cm x 10 cm at isocenter
- ▶ Tabletop material: Carbon-fiber  $\leq 1.4$  mm Al
- ▶ Mattress: standard  $\leq 0.6$  mm Al, thick  $\leq 1.0$  mm Al
- ▶ Different irradiation geometry:






# Interventional X-ray system

Table thickness crossed by the radiation is different for 0°, 30° and 90°.

$$d_2 > d_1 > d_3 = 0$$

| Tube Voltage (kV) | Additional filtration (mm Cu) | Ka (mGy) 0° | Ka (mGy) 30° | Ka (mGy) 90° |
|-------------------|-------------------------------|-------------|--------------|--------------|
| 70                | 0                             | 0.447       | 0.408        | 0.674        |
| 70                | 0.9                           | 0.021       | 0.020        | 0.029        |
| 125               | 0                             | 1.410       | 1.351        | 2.044        |
| 125               | 0.9                           | 0.299       | 0.286        | 0.390        |



# Conventional projection radiography

## ► Measurement geometry:





FIG. 3.2. Diagram of the measuring arrangement.



## ► Measurement equipment: XMM




# Conventional projection radiography

## ► Impact of field size:

| Tube voltage (kV) | Additional filtration (mm Cu) | HVL (mm Al) | Ka (mGy) 10x10 | Ka (mGy) 5x5 | Ka (mGy) 50x50 |
|-------------------|-------------------------------|-------------|----------------|--------------|----------------|
| 120               | 0                             | 4.91        | 1.1248         | 1.1036       | 1.2042         |
| 120               | 0.3                           | 8.88        | 0.4546         | 0.4473       | 0.4810         |

What is the reason? **Scatter radiation**

For example, scatter in the mirror of the X-ray tube



# Conventional projection radiography


► **Impact of the dosimeter position:**

| Tube voltage (kV) | Additional filtration (mm Cu) | HVL (mm Al) | Ka (mGy) 10x10 | Ka (mGy) 5x5 | Ka (mGy) 50x50 | Ka (mGy) 50x50 corner |
|-------------------|-------------------------------|-------------|----------------|--------------|----------------|-----------------------|
| 120               | 0                             | 4.91        | 1.1248         | 1.1036       | 1.2042         | 1.0937                |
| 120               | 0.3                           | 8.88        | 0.4546         | 0.4473       | 0.4810         | 0.3887                |

What is the reason? **Field inhomogeneity**

# Conventional projection radiography

## ► Impact of the dosimeter position:



In addition, the corner are not all the same. The values of air kerma measured in the corners are not equal, measurement in corner B is different from measurement in corner C.



# Thank you for your attention!



The project 22NRM01 TraMeXI has received funding from the European Partnership on Metrology, co-financed from the European Union's Horizon Europe Research and Innovation Programme and by the Participating States.  
Funded by the European Union. Views and opinions expressed are however those of the author(s) only and do not necessarily reflect those of the European Union or EURAMET. Neither the European Union nor the granting authority can be held responsible for them.

<http://tramexi.com>